

Thera Pearl Sports Pack



Reusable Hot & Cold Therapy

Rating: Not Rated Yet

Price

USD13.80

USD15.33

[Ask a question about this product](#)

Manufacturer [THERA°PEARL](#)

Description

Description:

The price of perfection is high, whether you pay at the gym or the spa. THERA°PEARL's Sports Pack rescues men and women from post-workout and post-treatment aches, pains, swelling and redness.

Compact and portable, the Doctor Designed Sports Pack is soft and pliable to adhere comfortably to your body. Drip and mess-free, it beats a melting ice pack or soggy hot towel any day! It's hypoallergenic, non-toxic and, best of all, reusable. For durability, cost and effectiveness, disposables can't match THERA°PEARL's value.

Chill it in the fridge or freezer for analgesic cold, or pop it in the microwave for penetrating heat. The pack holds its therapeutic temperature for a full 20 minutes so you get fast, soothing relief. It's perfect following any strenuous activity or accidental run-in with a table leg. Our Sports Pack is Perfect for:

- Sore Muscles/Arthritis
- Pre/Post Workout
- Sprains/Swelling & More!

Directions

THERA°PEARL Pack should be applied to skin for doctor recommended time of 20 minutes. Remove and wait at least 20 minutes. Reapply pack for 20 minutes, if needed.

HEAT THERAPY - Microwave for time indicated below and test pack temperature. Heat for additional 5 seconds as needed until desired temperature is reached.

- 700 Watt Microwave: 40 seconds
- 1000 Watt Microwave: 35 seconds
- 1250 Watt Microwave: 30 seconds

COLD THERAPY - Freeze for at least 2 hours.

Reviews

There are yet no reviews for this product.